

Габаритные размеры

ТИП	DNA	DNM	a	g	f	x1	x2	n1	n2	n3	h1	h2	h3	m1	z1	s	b	i1	i	L	d	t	u
100-150/2	150	100	289	296	470	629	589	435	355	275	250	280	280	70	70	20	80	166	135	110	48	51,8	14
100-150/3	150	100	393	296	470	733	693	435	355	275	250	280	280	70	70	20	80	166	135	110	48	51,8	14
100-150/4	150	100	497	296	470	837	797	435	355	275	250	280	280	70	70	20	80	166	135	110	48	51,8	14
100-150/5	150	100	601	296	470	941	901	435	355	275	250	280	280	70	70	20	80	166	135	110	48	51,8	14
100-150/6	150	100	705	296	470	1045	1005	435	355	275	250	280	280	70	70	20	80	166	135	110	48	51,8	14
100-150/7	150	100	809	296	470	1149	1109	435	355	275	250	280	280	70	70	20	80	166	135	110	48	51,8	14
100-150/8	150	100	913	296	470	1253	1213	435	355	275	250	280	280	70	70	20	80	166	135	110	48	51,8	14
100-150/9	150	100	1017	296	470	1357	1317	435	355	275	250	280	280	70	70	20	80	166	135	110	48	51,8	14
100-150/10	150	100	1121	296	470	1461	1421	435	355	275	250	280	280	70	70	20	80	166	135	110	48	51,8	14

